



SILLON DE DONACIÓN DE SANGRE - BLOOD DONATION ARMCHAIR



T4I

La línea de sillones de extracción Presvac, tiene todas las características técnicas necesarias para la extracción de sangre del donante, de manera confortable y segura, además de brindarle al operador una óptima y eficiente postura de trabajo.

- Accionamiento manual con mecanismos hidráulicos autofrenantes.
- Soporte lateral para agitador de extracción y elementos de donación.
- Totalmente articulado y de fácil limpieza.
- Dimensiones:
Sentado: ancho 80 cm - profundidad 103 cm - alto 136 cm.
Trendelenburg: ancho 80 cm - largo 174 cm - alto 79 cm.
- Estructura interior y cuerpo en acero.
- Ruedas giratorias con freno.

Presvac's blood donation armchairs are optimized with all the essential technical features to perform successful blood extractions in a safe and comfortable way, as well as providing the operator with an optimal and efficient working posture.

- Manual operation with hydraulic auto brake mechanisms.
- Lateral support to place extraction shaker and equipment.
- Completely articulated and easy to clean.
- Dimensions:
Sitting: width 80 cm - depth 103 cm - height 132 cm .
Trendelenburg: width 80 cm - length 174 cm - height 79 cm.
- Steel internal structure.
- Mounted on casters with braking system.

